

Anti-PP1R3A antibody (830-880) {BIOTIN} (STJ506060)

STJ506060

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PP1R3A (830-880) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5506
Gene Symbol	PPP1R3A
Uniprot ID	PPR3A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 830-880 on mouse Protein Phosphatase 1 Regulatory Subunit 3A.
Immunogen Region	830-880
Specificity	
Immunogen Sequence	