

Anti-SARS-CoV-2 membrane glycoprotein antibody {BIOTIN} (STJ506075)

STJ506075

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SARS-CoV-2 membrane glycoprotein is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken from unique amino acid sequence on membrane glycoprotein [Severe acute respiratory syndrome coronavirus 2].
Immunogen Region	
Specificity	
Immunogen Sequence	