

Anti-ABL2 antibody (1100-1150) (STJ506088)

STJ506088

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ABL2 (1100-1150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27
Gene Symbol	ABL2
Uniprot ID	ABL2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 1100-1150 on human Tyrosine-protein kinase ABL2/ARG.
Immunogen Region	1100-1150
Specificity	
Immunogen Sequence	