

**Anti-ABL2 antibody (1100-1150) {BIOTIN} (STJ506089)**

STJ506089

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-ABL2 (1100-1150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	BIOTIN
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	27
<b>Gene Symbol</b>	ABL2
<b>Uniprot ID</b>	ABL2_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 1100-1150 on human Tyrosine-protein kinase ABL2/ARG.
<b>Immunogen Region</b>	1100-1150
<b>Specificity</b>	
<b>Immunogen Sequence</b>	