

Anti-PTBP1 antibody (150-200) {BIOTIN} (STJ506092)

STJ506092

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PTBP1 (150-200) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5725
Gene Symbol	PTBP1
Uniprot ID	PTBP1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 150-200 on human Polypyrimidine tract-binding protein 1.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	