

Anti-MLKL antibody (400-450) {BIOTIN} (STJ506095)

STJ506095

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MLKL (400-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74568
Gene Symbol	Mkl
Uniprot ID	MLKL_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 400-450 on mouse Mixed lineage kinase domain-like protein MLKL.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	