

Anti-Phospho-PDE11A antibody (100-200) {BIOTIN} (STJ506101)

STJ506101

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PDE11A (100-200) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	
Formulation	
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	140928
Gene Symbol	Pde11a
Uniprot ID	PDE11_RAT
Immunogen	Synthetic peptide taken within amino acid region 100-200 on rat PDE11A protein. Phosphorylated Ser 163.
Immunogen Region	100-200
Specificity	
Immunogen Sequence	