

## Anti-Phospho-ARHGEF3 antibody (350-410) (STJ506109)

STJ506109

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-ARHGEF3 (350-410) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	50650
<b>Gene Symbol</b>	ARHGEF3
<b>Uniprot ID</b>	ARHG3_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 350-410 human Rho-Guanine nucleotide Exchange factor 3 protein.
<b>Immunogen</b>	350-410
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	