

Anti-Phospho-ARHGEF3 antibody (350-410) {BIOTIN} (STJ506110)

STJ506110

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-ARHGEF3 (350-410) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	50650
Gene Symbol	ARHGEF3
Uniprot ID	ARHG3_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 350-410 human Rho-Guanine nucleotide Exchange factor 3 protein.
Immunogen	350-410
Region	
Specificity	
Immunogen	
Sequence	