

Anti-CGRP antibody (350-410) (STJ506112)

STJ506112

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CGRP (350-410) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12310
Gene Symbol	Calca
Uniprot ID	CALCA_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 350-410 mouse Calcitonin gene-related peptide 1.
Immunogen Region	350-410
Specificity	
Immunogen Sequence	