

Anti-HYPE antibody (400-450) {BIOTIN} (STJ506116)

STJ506116

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HYPE (400-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11153
Gene Symbol	FICD
Uniprot ID	FICD_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 400-450 on mouse Protein adenyllyltransferase FICD.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	