

Anti-Profilin antibody (50-100) {FITC} (STJ506120)

STJ506120

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Profilin (50-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5216
Gene Symbol	PFN1
Uniprot ID	PROF1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 50-100 on human Profilin 1.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	