

## Anti-NHERF-1 antibody (255-305) (STJ506128)

STJ506128

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-NHERF-1 (255-305) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">9368</a>
<b>Gene Symbol</b>	<a href="#">SLC9A3R1</a>
<b>Uniprot ID</b>	<a href="#">NHRF1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 255-305 on human (+) /H (+) exchange regulatory cofactor NHERF1.
<b>Immunogen Region</b>	255-305
<b>Specificity</b>	
<b>Immunogen Sequence</b>	