

## Anti-TREML1 antibody (150-200) {BIOTIN} (STJ506132)

STJ506132

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-TREML1 (150-200) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	BIOTIN
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">340205</a>
<b>Gene Symbol</b>	<a href="#">TREML1</a>
<b>Uniprot ID</b>	<a href="#">TRML1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 150-200 on human Triggering Receptor Expressed on Myeloid cells-like 1.
<b>Immunogen Region</b>	150-200
<b>Specificity</b>	
<b>Immunogen Sequence</b>	