

Anti-TNK1 antibody (600-666) {BIOTIN} (STJ506150)

STJ506150

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TNK1 (600-666) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8711
Gene Symbol	TNK1
Uniprot ID	TNK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 600-666 on human Non-receptor tyrosine-protein kinase TNK1.
Immunogen Region	600-666
Specificity	
Immunogen Sequence	