

## Anti-TNK1 antibody (600-666) {FITC} (STJ506151)

STJ506151

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-TNK1 (600-666) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| <b>Applications</b>      | ELISA, IP, WB  |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human, Mouse, Rat  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.5-1.5 µg/µl   |
| <b>Conjugation</b>         | FITC  |
| <b>Purification</b>        | Affinity Purified   |
| <b>Dilution Range</b>      | WB: 1:500-1:5, 000<br>ELISA: 1:50, 000<br>IP: 1:200             |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | <a href="#">8711</a>   |
| <b>Gene Symbol</b>        | <a href="#">TNK1</a>   |
| <b>Uniprot ID</b>         | <a href="#">TNK1_HUMAN</a>   |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 600-666 on human Non-receptor tyrosine-protein kinase TNK1. |
| <b>Immunogen Region</b>   | 600-666  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |