

Anti-PARD6B antibody (340-372) {BIOTIN} (STJ506158)

STJ506158

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PARD6B (340-372) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84612
Gene Symbol	PARD6B
Uniprot ID	PAR6B_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 340-372 on human Partitioning Defective 6 homolog beta.
Immunogen Region	340-372
Specificity	
Immunogen Sequence	