

Anti-RHBDF2 antibody (1-50) {BIOTIN} (STJ506163)

STJ506163

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RHBDF2 (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79651
Gene Symbol	RHBDF2
Uniprot ID	RHBDF2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Inactive rhomboid protein 2.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	