

Anti-C1ORF115 antibody (75-125) {BIOTIN} (STJ506175)

STJ506175

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C1ORF115 (75-125) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79762
Gene Symbol	C1orf115
Uniprot ID	CA115_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 75-125 on human Chromosome 1 Open Reading Frame 115.
Immunogen Region	75-125
Specificity	
Immunogen Sequence	