

Anti-ERVW-1 antibody (480-538) {BIOTIN} (STJ506186)

STJ506186

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ERVW-1 (480-538) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	30816
Gene Symbol	ERVW-1
Uniprot ID	SYCY1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 480-538 on human ERVW-1 protein.
Immunogen Region	480-538
Specificity	
Immunogen Sequence	