

Anti-CXCR2 antibody (445-538) {BIOTIN} (STJ506204)

STJ506204

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CXCR2 (445-538) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12765
Gene Symbol	Cxcr2
Uniprot ID	CXCR2_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 445-538 on mouse C-X-C Chemokine Receptot type 2.
Immunogen Region	445-538
Specificity	
Immunogen Sequence	