

Anti-NEDD8 antibody (50-81) {BIOTIN} (STJ506210)

STJ506210

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NEDD8 (50-81) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4738
Gene Symbol	NEDD8
Uniprot ID	NEDD8_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 50-81 on human NEDD8 protein.
Immunogen Region	50-81
Specificity	
Immunogen Sequence	