

## Anti-IL18 antibody (37-193) {FITC} (STJ506220)

STJ506220

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-IL18 (37-193) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 Mu g/Mu l
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3606</a>
<b>Gene Symbol</b>	<a href="#">IL18</a>
<b>Uniprot ID</b>	<a href="#">IL18_HUMAN</a>
<b>Immunogen</b>	Recombinant protein corresponding to amino acids within region 37-193 on human IL18 protein.
<b>Immunogen Region</b>	37-193
<b>Specificity</b>	
<b>Immunogen Sequence</b>	