

Anti-SLC26A3 antibody (650-700) (STJ506233)

STJ506233

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-SLC26A3 (650-700) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5-1.5 Mu g/Mu l |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100 ICC: 1:50-1:100 IF: 1:50-1:100 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 1811 |
| Gene Symbol | SLC26A3 |
| Uniprot ID | S26A3_HUMAN |
| Immunogen | Synthetic peptide corresponding to amino acids within region 650-700 on human SCL26A3 protein. |
| Immunogen Region | 650-700 |
| Specificity | |
| Immunogen Sequence | |