

Anti-SLC26A6 antibody (700-759) (STJ506236)

STJ506236

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC26A6 (700-759) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100 ICC: 1:50-1:100 IF: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65010
Gene Symbol	SLC26A6
Uniprot ID	S26A6_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 700-759 on human SCL26A6 protein.
Immunogen Region	700-759
Specificity	
Immunogen Sequence	