

## Anti-SLC26A6 antibody (200-250) {BIOTIN} (STJ506237)

STJ506237

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-SLC26A6 (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 Mu g/Mu l
<b>Conjugation</b>	BIOTIN
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100 ICC: 1:50-1:100 IF: 1:50-1:100
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	65010
<b>Gene Symbol</b>	SLC26A6
<b>Uniprot ID</b>	S26A6_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 200-250 on rat 5-hydroxytryptamine receptor 1A protein.
<b>Immunogen Region</b>	200-250
<b>Specificity</b>	
<b>Immunogen Sequence</b>	