

Anti-EPS8 antibody (C-Term) (STJ70422)

STJ70422

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-EPS8 (C-Term) is suitable for use in ELISA, Western Blot and Immunohistochemistry research applications.
Applications	Pep-ELISA, WB, IHC
Host/Source	Goat
Reactivity	Human, Mouse, Rat, Dog, Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	IHC-Recommended concentration, 5-10µg/ml ELISA-antibody detection limit dilution 1:8000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2059
Gene Symbol	EPS8
Uniprot ID	EPS8_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	SGVESFDEGSSH



STJ70422 (5µg/ml) staining of paraffin embedded Human Breast Carcinoma. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.