

Anti-PALS2/VAM1/MPP6 antibody (N-Term) (STJ70484)

STJ70484

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-PALS2/VAM1/MPP6 (N-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB-Trf
Host/Source	Goat
Reactivity	Human, Mouse, Rat, Dog, Pig, Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	WB-0.3-1µg/ml IHC-5µg/ml ELISA-antibody detection limit dilution 1:128000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51678
Gene Symbol	PALS2
Uniprot ID	PALS2_HUMAN
Immunogen	
Immunogen Region	N-Term
Specificity	
Immunogen Sequence	QQVLENLTELPS



STJ70484 staining (0.3µg/ml) of 293 lysate (RIPA buffer, 35µg total protein per lane). Detected by western blot using chemiluminescence.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081