

Anti-FOXP4 antibody (N-Term) (STJ70568)

STJ70568

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-FOXP4 (N-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human, Mouse, Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	WB-1-3µg/ml ELISA-antibody detection limit dilution 1:4000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116113
Gene Symbol	FOXP4
Uniprot ID	FOXP4_HUMAN
Immunogen	
Immunogen Region	N-Term
Specificity	This antibody is expected to recognize all three reported isoforms (NP_001012426.1; NP_612466.1; NP_001012427.1).
Immunogen Sequence	VESASETIRSAAPSG