

Anti-CCL3 and CCL3L1 antibody (C-Term) (STJ70658)

STJ70658

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-CCL3 and CCL3L1 (C-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:32000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6348
Gene Symbol	CCL3
Uniprot ID	CCL3_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	Please note this product is expected to recognize the products of 2 different genes which are almost identical.
Immunogen Sequence	EWVQKYVSDLELSA