

Anti-Ghrelin preproprotein antibody (C-Term) (STJ70752)

STJ70752

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Ghrelin preproprotein (C-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	IHC-2-5µg/ml, 2ug/ml ELISA-antibody detection limit dilution 1:8000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51738
Gene Symbol	GHRL
Uniprot ID	GHRL_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	This antibody is expected to recognise all reported isoforms.
Immunogen Sequence	QDILWEEAKEAPADK