

## Anti-Foxi2-Mouse antibody (C-Term) (STJ70812)

STJ70812

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Goat polyclonal antibody anti-Foxi2-Mouse (C-Term) is suitable for use in ELISA research applications. |
| <b>Applications</b>      | Pep-ELISA  |
| <b>Host/Source</b>       | Goat   |
| <b>Reactivity</b>        | Mouse, Rat   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.5 mg/mL   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| <b>Dilution Range</b>      | ELISA-antibody detection limit dilution 1:8000.   |
| <b>Formulation</b>         | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20 on receipt and minimise freeze-thaw cycles.  |

### TARGET INFORMATION

|                           |             |
|---------------------------|-------------|
| <b>Gene ID</b>            |             |
| <b>Gene Symbol</b>        |             |
| <b>Uniprot ID</b>         |             |
| <b>Immunogen</b>          |             |
| <b>Immunogen Region</b>   | C-Term      |
| <b>Specificity</b>        |             |
| <b>Immunogen Sequence</b> | GHLIYSRDGTE |