

Anti-ETV1 antibody (Internal) (STJ70847)

STJ70847

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-ETV1 (Internal) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human, Mouse, Rat, Dog

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:8000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2115
Gene Symbol	ETV1
Uniprot ID	ETV1_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	This antibody is expected to recognise all reported isoforms. Variants (NP_001156620.1; NP_001156621.1) encode the same isoform.
Immunogen Sequence	PLKIKKKEPHSPCSE