

Anti-Rabring 7/BCA2 antibody (Internal) (STJ70879)

STJ70879

GENERAL INFORMATION


| | |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Goat polyclonal antibody anti-Rabring 7/BCA2 (Internal) is suitable for use in ELISA, Western Blot, Immunohistochemistry and Flow Cytometry research applications. |
| Applications | Pep-ELISA, WB, IHC, FC |
| Host/Source | Goat |
| Reactivity | Human, Mouse, Rat, Dog |

PRODUCT PROPERTIES

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Dilution Range | WB-0.3-1µg/ml ELISA-antibody detection limit dilution 1:32000. |
| Formulation | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Isotype | IgG |
| Storage Instruction | Store at -20 on receipt and minimise freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 27246 |
| Gene Symbol | RNF115 |
| Uniprot ID | RN115_HUMAN |
| Immunogen | |
| Immunogen Region | Internal |
| Specificity | |
| Immunogen Sequence | DQDNRANERGHQTH |



STJ70879 (0.3µg/ml) staining of HEPG2 cell lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.