

Anti-X-arrestin/arrestin 3 antibody (C-Term) (STJ71085)

STJ71085

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-X-arrestin/arrestin 3 (C-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	IHC-2.5µg/ml ELISA-antibody detection limit dilution 1:16000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	407
Gene Symbol	ARR3
Uniprot ID	ARRC_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	RKGEEESQKAVE



STJ71085 (2.5µg/ml) staining of paraffin embedded Human Retina. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.

STJ71085 Negative Control showing staining of paraffin embedded Human Retina, with no primary antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081