

Anti-PPP2R3A antibody (721-31) (STJ71329)

STJ71329

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-PPP2R3A (721-31) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:1000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5523
Gene Symbol	PPP2R3A
Uniprot ID	P2R3A_HUMAN
Immunogen	
Immunogen Region	721-31
Specificity	This antibody is expected to recognise both reported isoforms (NP_002709.2; NP_871626.1)
Immunogen Sequence	TCSNHEQTLRSR