

Anti-Urocortin 3 antibody (Internal) (STJ71448)

STJ71448

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Urocortin 3 (Internal) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:16000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114131
Gene Symbol	UCN3
Uniprot ID	UCN3_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	KPRQDTAKSPHRTK



HEK293 lysate (10µg protein in RIPA buffer) over expressing Human UCN3 with DYKDDDDK tag probed with STJ71448 (1µg/ml) in Lane A and probed with anti-DYKDDDDK tag (1/3000) in lane C. Mock-transfected HEK293 probed with STJ71448 (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081