

Anti-MUPCDH antibody (C-Term) (STJ71484)

STJ71484

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-MUPCDH (C-Term) is suitable for use in ELISA, Western Blot and Flow Cytometry research applications.
Applications	Pep-ELISA, WB, FC
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:8000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53841
Gene Symbol	CDHR5
Uniprot ID	CDHR5_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	This antibody is expected to recognise both reported isoforms (NP_068743.2; NP_112554.2).
Immunogen Sequence	PYDAPGGDDSY