

## Anti-TUSC2/FUS1 antibody (Internal) (STJ71528)

STJ71528

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Goat polyclonal antibody anti-TUSC2/FUS1 (Internal) is suitable for use in ELISA and Western Blot research applications. |
| <b>Applications</b>      | Pep-ELISA, WB  |
| <b>Host/Source</b>       | Goat   |
| <b>Reactivity</b>        | Human, Mouse, Rat  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.5 mg/mL   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| <b>Dilution Range</b>      | ELISA-antibody detection limit dilution 1:2000.   |
| <b>Formulation</b>         | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20 on receipt and minimise freeze-thaw cycles.  |

### TARGET INFORMATION

|                           |                |
|---------------------------|----------------|
| <b>Gene ID</b>            | 11334          |
| <b>Gene Symbol</b>        | TUSC2          |
| <b>Uniprot ID</b>         | TUSC2_HUMAN    |
| <b>Immunogen</b>          |                |
| <b>Immunogen Region</b>   | Internal       |
| <b>Specificity</b>        |                |
| <b>Immunogen Sequence</b> | TKNGQKRAKLRRVH |