

Anti-CHRFAM7A/CHRNA7-FAM7A antibody (N-Term) (STJ71540)

STJ71540

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-CHRFAM7A/CHRNA7-FAM7A (N-Term) is suitable for use in ELISA and Immunohistochemistry research applications.
Applications	Pep-ELISA, IHC
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:1000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage	Store at -20 on receipt and minimise freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	89832
Gene Symbol	CHRFAM7A
Uniprot ID	CRFM7_HUMAN
Immunogen	
Immunogen	N-Term
Region	
Specificity	This antibody is expected to recognize reported isoform NP_647536.1 only.
Immunogen	QKYCIYQHFQFQ
Sequence	