

Anti-Neurotrophin 3 antibody (Internal) (STJ71567)

STJ71567

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Neurotrophin 3 (Internal) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human, Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:1000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4908
Gene Symbol	NTF3
Uniprot ID	NTF3_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	This antibody is expected to recognise both isoforms (NP_001096124.1; NP_002518.1)
Immunogen Sequence	TSRRKRYAEHKSH