

Anti-TFAP2D antibody (Internal) (STJ71712)

STJ71712

GENERAL INFORMATION

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Goat polyclonal antibody anti-TFAP2D (Internal) is suitable for use in ELISA and Western Blot research applications. |
| Applications | Pep-ELISA, WB |
| Host/Source | Goat |
| Reactivity | Human, Mouse, Rat, Horse, Dog, Cow |

PRODUCT PROPERTIES

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Dilution Range | ELISA-antibody detection limit dilution 1:64000. |
| Formulation | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Isotype | IgG |
| Storage Instruction | Store at -20 on receipt and minimise freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------|
| Gene ID | 83741 |
| Gene Symbol | TFAP2D |
| Uniprot ID | AP2D_HUMAN |
| Immunogen | |
| Immunogen Region | Internal |
| Specificity | |
| Immunogen Sequence | HHQSFHYEFQHSHP |