

Anti-BAT3/BAG6 antibody (Internal) (STJ71777)

STJ71777

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-BAT3/BAG6 (Internal) is suitable for use in Western Blot research applications.
Applications	WB-Trf
Host/Source	Goat
Reactivity	Human, Mouse, Rat, Dog, Cow, Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	IHC-Recommended concentration, 3-5µg/ml ELISA-antibody detection limit dilution 1:128000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7917
Gene Symbol	BAG6
Uniprot ID	BAG6_HUMAN
Immunogen Region	Internal
Specificity	This antibody is expected to recognize both reported isoforms (NP_004630.3; NP_542433.1). Reported variants representing identical protein: NP_542434.1, NP_001092004.1, NP_542433.1 (isoform b).
Immunogen Sequence	PETNAPNHPSPAEY