

Anti-TCEAL1 antibody (N-Term) (STJ71878)

STJ71878

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-TCEAL1 (N-Term) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:8000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9338
Gene Symbol	TCEAL1
Uniprot ID	TCAL1_HUMAN
Immunogen	
Immunogen Region	N-Term
Specificity	Reported variants represent identical protein: NP_001006641.1, NP_001006640.1, NP_004771.2.
Immunogen Sequence	DKPRKENEEEPQ