

Anti-MOZ/KAT6A antibody (Internal) (STJ71939)

STJ71939

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-MOZ/KAT6A (Internal) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	IHC-3.75µg/ml ELISA-antibody detection limit dilution 1:16000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7994
Gene Symbol	KAT6A
Uniprot ID	KAT6A_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	Reported variants NP_001092882.1, NP_006757.2, NP_001092883.1 represent identical protein.
Immunogen Sequence	EKIKDKKEETELD