

Anti-C7orf50 antibody (Internal) (STJ72079)

STJ72079

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-C7orf50 (Internal) is suitable for use in ELISA, Western Blot, Immunohistochemistry and Immunocytochemistry research applications.
Applications	Pep-ELISA, WB, IHC, ICC
Host/Source	Goat
Reactivity	Human, Dog

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:32000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage	Store at -20 on receipt and minimise freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84310
Gene Symbol	C7orf50
Uniprot ID	CG050_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	Reported variants represent identical protein: NP_001127867.1, NP_001127868.1, NP_115726.1
Immunogen Sequence	EGLQGRARELTVQK