

Anti-FOXI3/Forkhead box I3 antibody (Internal) (STJ72091)

STJ72091

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-FOXI3/Forkhead box I3 (Internal) is suitable for use in ELISA, Immunohistochemistry and Immunocytochemistry research applications.
Applications	Pep-ELISA, IHC, ICC
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	IHC-5-10µg/ml ELISA-antibody detection limit dilution 1:4000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	344167
Gene Symbol	FOXI3
Uniprot ID	FOXI3_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	SPTSISEASADTLQ