

Anti-UBP1 antibody (Internal) (STJ72115)

STJ72115

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Goat polyclonal antibody anti-UBP1 (Internal) is suitable for use in ELISA research applications. |
| Applications | Pep-ELISA |
| Host/Source | Goat |
| Reactivity | Human, Mouse, Rat, Dog |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Dilution Range | IHC-10µg/ml ELISA-antibody detection limit dilution 1:2000. |
| Formulation | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA |
| Isotype | IgG |
| Storage Instruction | Store at-20 on receipt and minimise freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 7342 |
| Gene Symbol | UBP1 |
| Uniprot ID | UBIP1_HUMAN |
| Immunogen | |
| Immunogen Region | Internal |
| Specificity | Reported variants represent identical protein: NP_001121633.1, NP_055332.3. |
| Immunogen Sequence | PNHQGDGASQTSGE |