

Anti-FANCE antibody (Internal) (STJ72281)

STJ72281

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-FANCE (Internal) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB-Trf
Host/Source	Goat
Reactivity	Human, Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	IHC-Recommended concentration, 3-6µg/ml ELISA-antibody detection limit dilution 1:32000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at-20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2178
Gene Symbol	FANCE
Uniprot ID	FANCE_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	EHSLESLADGGS