

## Anti-Desmoplakin antibody (1546-1558) (STJ72284)

STJ72284

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Goat polyclonal antibody anti-Desmoplakin (1546-1558) is suitable for use in ELISA, Western Blot and Immunohistochemistry research applications.
<b>Applications</b>	Pep-ELISA, WB, IHC
<b>Host/Source</b>	Goat
<b>Reactivity</b>	Human, Mouse, Rat, Cow

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Dilution Range</b>	ELISA-antibody detection limit dilution 1:4000.
<b>Formulation</b>	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20 on receipt and minimise freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">1832</a>
<b>Gene Symbol</b>	<a href="#">DSP</a>
<b>Uniprot ID</b>	<a href="#">DESP_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	1546-1558
<b>Region</b>	
<b>Specificity</b>	This antibody is expected to recognize reported isoform I (NP_004406.2) only.
<b>Immunogen</b>	SQERTVKDQDITR
<b>Sequence</b>	